

Honeywell

CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
lab.honeywell.com

Brand: Honeywell Riedel-de Haën

Product: 34850

Lot No: O173M

Production date: 21.Jun.2024

Expiry Date: 08.Dec.2026

Acetone CHROMASOLV, for HPLC, 99.8%

Parameter	Specification	Units	Result	Units
assay (GC)	min. 99.8	%	99.98	%
non-volatile matter	max. 0.0005	%	<0.0005	%
water (Karl Fischer)	max. 0.1	%	0.03	%
free acid (as CH3COOH)	max. 0.002	%	<0.002	%
transmittance at 330 nm	min. 15	%	>15	%
transmittance at 335 nm	min. 50	%	>50	%
transmittance at 340 nm	min. 80	%	>80	%
transmittance from 350 nm	min. 98	%	>98	%



filtered via 0,2µm pore size filter

QC Release Date: 15.Aug.2024

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell

Quality Management

This document was produced electronically and is valid without a signature