



CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
lab.honeywell.com

Brand: Honeywell Riedel-de Haën
Product: 34854
Lot No: O186M
Production date: 04.Jul.2024
Expiry Date: 24.Jun.2026

Chloroform CHROMASOLV, for HPLC, 99.8%, amylene stabilized

Parameter	Specification	Units	Result	Units
assay (GC)	min. 99.80	%	99.91	%
non-volatile matter	max. 0.0005	%	<0.0005	%
water (Karl Fischer)	max. 0.01	%	<0.01	%
free acid (as HCl)	max. 0.001	%	<0.001	%
absorbance at 250 nm	max. 1.00		<1.00	
absorbance at 260 nm	max. 0.30		<0.30	
absorbance 280 nm	max. 0.005		<0.005	
IR spectrum	complying		complying	



filtered via 0,2µm pore size filter

QC Release Date: 31.Jul.2024

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell

Quality Management

This document was produced electronically and is valid without a signature